

Cochlear™ Nucleus® System Accessories

Maximize Listening



Mini Microphone 2+

Wirelessly hear a speaker from a distance or conversations in noisy listening situations such as in meetings, in the car, at restaurants, or in a classroom.

N6 N7 K

- Clip or place next to a sound source, and it will transmit the sound wirelessly over a distance up to 82 feet
- Works with one or two Nucleus 7, Nucleus 6 or Kanso® Sound Processors
- Connects to FM systems
- Plug into any other sound device like an MP3 Player and more to bring the sound directly to you
- Wireless alternative to the Lapel Microphone for the Nucleus 6 Sound Processor



Euro Accessory Adaptor

Connect your FM system directly to your Nucleus 6 Sound Processor and the voice of a speaker transmits to your sound processor wirelessly in the classroom or in noisy listening situations.

N6

- Some ear-level FM receivers such as the Phonak MLxS cannot connect directly to your Nucleus 6 Sound Processor without the Euro Accessory Adaptor
- The Euro Accessory Adaptor provides a convenient interface to connect these ear-level FM receivers directly to your sound processor
- Use with a third-party FM System



Lapel Microphone

Hear conversations more easily in noisy situations such as in a restaurant or while traveling in an automobile.

N6

- Clip the microphone close to whomever you're chatting with to capture sounds coming from various directions
- Use with the Nucleus 6 Sound Processor with an accessory port



Loop Booster

Minimize interference by placing the loop booster away from the processor when using your Nucleus 6 Sound Processor with accessory port.

N6

Cochlear™ Nucleus® System Accessories

Enjoy Sports and Activities



Cochlear Aqua+ Kit for Nucleus 6



Cochlear Aqua+ for Nucleus 7

Swim underwater, in salt water, and completely submerged without the need for any cumbersome cables, cases, or a special sound processor specifically for water activities***

N6 N7 K

- Ideal for those who wish to have additional water resistance for strenuous activities, too
- Can be submerged in up to 10 feet for 2 hours
- Can be reused up to 50 times
- Use the Aqua+ Compact for Nucleus 7 with compact rechargeable battery modules



Cochlear Aqua+ for Kanso



Cochlear Aqua+ Compact for Nucleus 7



Aqua Accessory

Enjoy water activities with an economical and disposable water accessory option.***

N6

- Ideal for those who love to swim underwater or wish to have additional water resistance for strenuous activities
- Insert the sound processor in the Aqua Accessory plastic enclosure with zip lock seal
- Alternative to the Nucleus Aqua+

* The Nucleus Aqua+ accessory is IP68 rated and is approved for use with rechargeable batteries and is not approved for use with the acoustic component. This water protection rating means that the sound processor with the Aqua+ can be continuously submerged under water to a depth of 3 meters (9 feet and 9 inches) for up to 2 hours.

** The Kanso Sound Processor with the Aqua+ is water resistant to level IP68 of the International Standard IEC60529. This water protection rating means that the sound processor with the Aqua+ can be continuously submerged under water to a depth of 3 meters (9 feet and 9 inches) for up to 2 hours. This water protection only applies when you use the Aqua+ and LR44 alkaline or nickel metal hydride disposable batteries.

*** The Aqua Accessory has achieved quality waterproof ratings in depths up to 4 meters for duration of 2 hours.



Hear now. And always

Cochlear®

Cochlear™ Nucleus® System Accessories

Listen to Music, TV and Other Devices



Direct Streaming

Use your compatible Apple® or Android™ device to stream sound directly to your Nucleus 7 Sound Processor without the need for any attachments to the sound processor.†

N7

- You can enjoy your favorite movies, TV shows, apps, music, get GPS directions—even a FaceTime® call with audio streamed directly to your sound processor
- Designed to provide a better experience and fit seamlessly into your lifestyle



Phone Clip

Wirelessly stream hands-free phone calls and audio directly to your sound processor from Bluetooth® enabled devices.

N6 N7 K

- Works with one or two Nucleus 7, Nucleus 6 or Kanso Sound Processors
- Answer and adjust the volume just by pushing a button
- Transmits audio wirelessly over a distance up to 23 feet
- Wireless alternative option to the Unilateral and Bilateral Personal Audio Cables and the Portable Phone Cable for the Nucleus 6/7 Sound Processors



TV Streamer

Enjoy sound directly from your TV, stereo or other audio sources wirelessly without turning up the volume and disturbing others.

N6 N7 K

- Works with one or two Nucleus 7, Nucleus 6 or Kanso Sound Processors
- Move freely, up to 23 feet away, while listening to the TV
- Wireless alternative option to the Mains Isolation Cable for the Nucleus 6/7 Sound Processors



Unilateral Personal Audio Cable

Hear music and other audio from electronic gadgets with a clearer, direct-to-sound processor connection when you use one Nucleus 6 Sound Processor.

N6

- Connect a battery-powered sound source (e.g., an iPod®) with a 3.5mm headphone jack to your one sound processor
- Produces monophonic sound
- 120cm/47" length
- Use with the Nucleus 6 Sound Processor with an accessory port



Freedom® Accessory Adaptor

Enjoy direct-to-sound processor connection with your previous sound processor generation audio cables.

N6

- Plug the Freedom Accessory adaptor into the Nucleus 6 Sound Processor
- Use with the Nucleus 6 Sound Processor with an accessory port
- Plug Nucleus Freedom Audio Cables into the Freedom Adaptor



Bilateral Personal Audio Cable

Hear music and other audio from electronic gadgets with a direct-to-sound processor connection for a stereo hearing experience when you use two Nucleus 6 Sound Processors.

N6

- Connect a battery-powered sound source with a 3.5mm jack to your two Nucleus 6 Sound Processors using just one cable – for surround sound enjoyment
- Produces either monophonic or stereophonic sound, depending on your sound system
- Use with the Nucleus 6 Sound Processor with an accessory port
- 120cm/47" length



Mains Isolation Cable

Enjoy hearing the TV, computer, and other electronic devices with a direct-to-sound processor connection.

N6

- Connect your sound processor to a "mains-powered" sound source with a 3.5mm headphone jack that plugs into a power outlet (e.g. TV, computer, etc.)
- 3m/9'10" length
- Incorporates surge protection
- and will help protect your Nucleus 6 Sound Processor
- **Important:** This cable is used with your Personal Audio Cable or Bilateral Audio Cable
- If your Nucleus 6 Sound Processor is equipped with an accessory port



Portable Phone Cable

Hear music and other audio from electronic gadgets with a clearer, direct-to-sound processor connection when you use one Nucleus 6 Sound Processor.

N6

- Connect a battery-powered sound source with a 2.5mm jack (e.g. some types of mobile phones)
- Does not have a microphone, but uses the phone's built-in microphone (i.e. the piece you talk into)
- Produces monophonic sound
- 120cm/47" length
- Use with the Nucleus 6 Sound Processor with an accessory port

Customize/Wear



LiteWear™

LiteWear provides more comfort and less weight on your ear by having the battery module worn on your body.

N6

- Great for children or anyone interested in reducing the sound processor weight on their ears
- Gives you the option of connecting the battery module to a piece of clothing



Mic Locks

Hold your sound processor more securely and comfortably on the ear especially when using the Nucleus Aqua+.

N6

- A clear band that encloses around the base of the sound processor, wraps around the ear, and connects to the earhook
- Easily customized to fit any size ear to provide a more secure fit, tamper resistance, and comfort
- Works with the Cochlear Nucleus Aqua+




Sound Processor and Coil Covers

Protect and style your sound processor and coil with 12-personality plus covers to make it your own.

N6

- The sound processor and coil covers attach in a quick snap
- Available in standard and compact sizes, 12 fun colors and pattern



Personalization Stickers

Show your personal style with a selection of stickers you can apply directly to your sound processor.

N6 N7 K

- A single sticker sheet includes a fun assortment of 50 sticker designs including koala, pizza and diamond.



Microphone Protectors

Enjoy crisper hearing with fresh microphone protectors designed to optimize your hearing performance.

N6

- It's important to have microphone protectors in place at all times
- Change both microphone protectors once every three months (*or more frequently if you live in a humid environment or notice a difference in your sound quality*)



Snugfit™

Easily tighten the sound processor on you or your child's ear for added flexibility, security and peace of mind.

N6

- Ideal for those who engage in physical activities where the sound processor is at risk of falling off the ear
- Available in Small, Medium, and Large; 4 color options
- May help those with small ears or who need an extra level of security or comfort than an earhook alone can provide



Hugfit™

Hugfit™ provides an additional retention option to help you wear your Nucleus 7 Sound Processor more securely.

N7

- For extra peace of mind, the Cochlear Hugfit holds a child's sound processor in place more securely than an earhook alone
- Available in sizes XXSmall, XSmall and Small



Microphone Covers

You'll enjoy crisper hearing with a set of new Nucleus 7 Microphone Covers, designed to optimize your hearing performance.

N7

- You should replace your Cochlear Nucleus 7 processor's microphone covers every three months, or sooner if you live in a high humidity climate
- The Microphone protectors are built into the Nucleus 7 Microphone cover, which is completely replaced
- **Important:** Keep your microphone protectors in place at all times



Koala Clip

For extra peace of mind, the Koala Clip is an additional retention device for children.

N7

- Children using the Nucleus 7 may benefit from using one of Cochlear's loss prevention accessories.**
- The Cochlear Koala Clip allows for the Nucleus 7 Sound Processor to be attached to a child's clothing until they can wear their processor on their ear



Microphone Protectors

Enjoy crisper hearing with microphone covers designed to optimize our hearing performance.

K

- You should replace your Kanso Microphone Protector every 3 months, or sooner if you live in a humid or dusty climate.
- **Important:** Keep your microphone protectors in place at all times.

** Retention lines longer than the Safety Line (standard length) are not recommended for use by children as they may present a risk of strangulation.

Customize/Wear



Nucleus 7 Pediatric Headband

Hold the Nucleus 7 Sound Processor comfortably and securely in place for added peace of mind.

N7

- Perfect for securing the Nucleus 7 Sound Processor during physical activities where the sound processor is at risk of falling off
- Available in XSmall and Small
- Available in the color black



Cochlear Headband

Hold the Kanso Sound Processor comfortably and securely in place for added peace of mind.

K

- Perfect for securing the Kanso Sound Processor during physical activities where the sound processor is at risk of falling off
- Available in XXSmall, XSmall, Small, Medium and Large
- Available in 5 colors



Long Safety Line with Alligator Clip

Easily secure the Kanso Sound Processor to an article of clothing to prevent loss.

K

- A clear safety line connects the sound processor to an article of clothing with a snap-lock clip for added retention
- Long safety line available – 26 cm (approx. 10.2")
- Works with the Aqua+ for Kanso
- Standard safety line available – 21 cm (approx. 8.3")



Short Safety Line with Hair Clip

For a more discreet option secure the Kanso Sound Processor with the short safety line and hair clip.

K

- A clear safety line connects the sound processor to the hair clip for a discreet retention option
- Short safety line available – 8 cm (approx. 3.1")
- Hair clip available in 3 colors
- Works with the Aqua+ for Kanso
- Includes Pack of 5 Hair Clips

Charge Your Battery



Cochlear Y Battery Charger

Charge a standard or compact rechargeable battery module with these Nucleus 7 battery chargers. Both can either be powered directly from a USB port (e.g. computer) or, when used with the Cochlear USB Power Adaptor, by a wall power outlet.

N7

- The Cochlear Nucleus 7 Y Battery Charger lets you charge two battery modules at a time
- Both chargers come with a USB cover and the Cochlear USB Power Adapter
- The Cochlear Nucleus 7 USB Battery Charger can charge a single battery module and is small enough to keep with you while you're on the go



Cochlear USB Battery Charger



Nucleus 6 Sound Processor Battery Charger

The Nucleus 6 Sound Processor Battery Charger and Battery Charging Kit allow you to charge your rechargeable batteries.

N6

- This charger base comes with four sockets, so you can charge up to four battery modules at the same time
- The Nucleus Battery Charging Kit consists of the wall adapter and plugs required to use the Battery Charger base
- An empty compact or standard battery module can be fully recharged in about four hours



Nucleus Battery Charging Kit

Monitor Your Hearing



CR310 Remote Control

Control and monitor True Wireless™ devices using the CR310 Remote Control.

N7

- Once paired with your Nucleus 7 Sound Processor, the Mini Microphone and TV Streamer can be controlled directly from your processor, using the CR310 Remote Control or with the Nucleus Smart App



Nucleus Smart App

Use the Nucleus Smart App† to control and adjust sound processor settings, track hearing information and battery levels from your compatible Apple® or Android™ device.

N7

- Monitor the status of your sound processor, adjust the mixing ratio or volume, and use the Hearing Tracker all from your compatible device without carrying additional equipment
- Find a misplaced or lost Nucleus 7 Sound Processor with the 'Find my Processor' function



Monitor Earphones

Use the monitor earphones to check the Nucleus 7 Sound Processor of a recipient in your care.

N7

- Caregivers can use the Monitor Earphone Adaptor to check that a Cochlear implant recipient is receiving sound, and that functions such as telecoil or wireless accessories are working properly
- The Monitor Earphone Adaptor fits between the rechargeable battery module and your processing unit



Monitor Earphone Adaptor

Provides the ability for you to hear whether your child's or loved one's sound processor is working.

N6

- Connect the monitor earphones to the sound processor and place the earbuds in the caregiver's ears to check that the sound processor is working
- Also, use the accessory socket built into the monitor earphones to check if wired accessories such as the Personal Audio Cable and Euro Accessory Adaptor works
- 160cm/63" length



CR230 Remote Assistant

Put hearing control in your own hands and adjust settings at your fingertips with this advanced Remote Assistant.

N6

K

- A great tool for parents and caregivers who wish to monitor their loved one's hearing or for those who like to have additional control over their sound processor
- Once wirelessly paired, you can troubleshoot, view the sound processor's status, adjust volume/sensitivity, switch programs, and use telecoil
- Uses a built-in rechargeable battery



CR210 Remote Control

Small-sized remote control helps you change volume and control your hearing at the push of a button discreetly and simply.

N6

K

- Attach it to your keychain or even put it in your pocket
- Once wirelessly paired, you can change volume/sensitivity, switch programs, or use telecoil
- Offers longer battery life than the CR230 Remote Assistant and uses disposable batteries




www.Cochlear.com/US

Cochlear Americas

13059 East Peakview Avenue
Centennial, CO 80111 USA
Telephone: 1 303 790 9010
Support: 1 800 483 3123

Cochlear Canada Inc.

2500-120 Adelaide Street West
Toronto, ON M5H 1T1 Canada
Telephone: 1 800 483 3123

-  Compatible with the Nucleus® 7 Sound Processor
-  Compatible with the Nucleus® 6 Sound Processor
-  Compatible with the Nucleus® Kanso® Sound Processor

©Cochlear Limited 2019. All rights reserved. Hear now. And always and other trademarks and registered trademarks are the property of Cochlear Limited or Cochlear Bone Anchored Solutions AB. The names of actual companies and products mentioned herein may be the trademarks of their respective owners.

Apple, the Apple logo, FaceTime, Made for iPad logo, Made for iPhone logo, Made for iPod logo, iPhone, iPad Pro, iPad Air, iPad mini, iPad and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries.

Android and Google Play are registered trademarks of Google Inc.

The Cochlear Nucleus Smart App is available on App Store and Google Play. For complete smartphone compatibility information, please visit www.cochlear.com/compatibility.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Cochlear is under license.



FUN2207 ISS6 OCT19

Hear now. And always

