

Previous Processors Compared to the

Cochlear™ Baha® 5 Sound Processor



Sound Processors	Baha® Divino®	Baha® 3 System	Baha® 4 System	Baha® 5 System
Wearing Style	Headworn	Headworn	Headworn	Headworn
Year	2004	2009	2013	2015
Size (L x W x H)	27 x 20 x 11mm	30 x 21 x 12mm	30 x 21 x 12mm	26 x 19 x 12mm
Weight	11g	11.3g (excl. battery)	11.6g (excl. battery)	9.8g (excl. battery)
Processor Color Options	4	6	5	5
Processing Features				
SmartSound® iQ	No	No	No	Yes
Bluetooth® Smart Technology	No	No	No	Yes
# of Channels	1	12	17	17
Programs	2	3	4	34*
Signal Processing	Digital	Digital	Wireless/Digital	Wireless/Digital
Wireless Streaming	No	No	Yes	Yes
Feedback Processing	No	Cancellation	Dual-Track Management	Dual-Track Management
Datalogging	No	Yes	Yes	Yes
2.4 GHz Wireless	No	No	Yes	Yes
Made for iPhone®	No	No	No	Yes
Smart App	No	No	No	Yes
Control App	No	No	Yes	Yes
Geotagging	No	No	No	Yes
Mixing Ratio Control	No	Yes†	Yes	Yes
"Locate Device" Function	No	No	No	Yes
Low Battery Warning	No	Yes	No	Yes
Treble and Bass Adjust	No	No	No	Yes

Hear now. And always





Sound Processors	Baha® Divino®	Baha® 3 System	Baha® 4 System	Baha® 5 System
Accessories				
Tamper-proof Battery Doors	No	Yes	Yes	Yes
Sticker Accessories	Yes	Yes	Yes	Yes
Direct Audio Input	Yes	Yes	Yes	No
Battery Door Color Options	No	12	12	1
Wireless Remote Control	No	No	Yes	Yes
Wireless Bluetooth-Enabled Phone Clip	No	No	Yes	Yes
Wireless Mini Mic	No	No	Yes	Yes
Wireless TV Streamer	No	No	Yes	Yes
Telecoil	Available	Available	Available	No
FM	Yes	Yes	Yes	Available [†]
Dust and Moisture Protection	None	Multi-layered microphone cover linings	Multi-layered microphone cover linings	Multi-layered microphone cover linings
Active Sound Programs Features				
Wind	No	Yes	Yes	Yes
Music	No	Yes	Yes	Yes
Noise	No	Yes	Yes	Yes
Sound Environment Management	Manual	Manual	Automatic/Adaptive AND Manual	Automatic/Adaptive AND Manual

^{*}Using the Smart App's Geotagging function

[†]With the use of a body-worn receiver and Mini Microphone

[‡]With DAI and microphone input

The Cochlear Baha 5 Sound Processor is compatible with iPhone 6 Plus, iPhone 6, iPhone 5s, iPhone 5c, iPhone 5, iPad Air 2, iPad Air, iPad (4th generation), iPad mini 3, iPad mini 2, iPad mini and iPod touch (5th generation) using iOS 7* or later.

The Cochlear Baha 5 Smart App is compatible with iPhone 6 Plus, iPhone 6, iPhone 5s, iPhone 5c, iPhone 5, iPad Air 2, iPad Air, iPad (4th generation), iPad mini 3, iPad mini 2, iPad mini and iPod touch (5th generation) using iOS 8.1 or later.

www.Cochlear.com/US

Cochlear Americas
13059 East Peakview Avenue
Centennial, CO 80111 USA
Telephone: 1 303 790 9010
Support: 1 800 483 3123

Cochlear Canada Inc.
2500-120 Adelaide Street West
Toronto, ON M5H 1T1 Canada
Telephone: 1 416 972 5082
Fax: 1 416 972 5083

©Cochlear Limited 2015. All rights reserved. Hear now. And always and other trademarks and registered trademarks are the property of Cochlear Limited or Cochlear Bone Anchored Solutions AB. The names of actual companies and products mentioned herein may be the trademarks of their respective owners.



BUN267 ISS3 DEC15

