

Surgical Instruments/Components for Surgery

Reusable Instruments Requiring Sterilization

The following instruments need to be sterilized. Two sets are needed if more than one surgery is planned, so that the second surgery can begin as soon as the first one is finished. Then, the first set is re-sterilized during the second surgery so it can be used for the third surgery.

New Instruments for Baha Attract			
93571	Indicator for Baha Attract	Non-sterile	Reusable
94071	Implant magnet template	Non-sterile	Reusable
93572	Bone bed indicator	Non-sterile	Reusable
93570	Soft tissue gauge, 6mm	Non-sterile	Reusable

Note:

- Thin hypodermic needles and a clamp (mosquito) are needed for measuring of soft tissue thickness prior to incision.
- A normal high speed otological drill unit and diamond burs are needed as back-up in case bone polishing is needed.

Existing Baha Instruments Needed for Baha Attract			
92146	Baha instrument cassette	Non-sterile	Reusable
92142	Implant inserter	Non-sterile	Reusable
90381	Machine screwdriver Unigrip, 25mm	Non-sterile	Reusable
90469	Screwdriver Unigrip, 95mm	Non-sterile	Reusable
90944	Raspatorium	Non-sterile	Reusable
90943	Dissector	Non-sterile	Reusable
92143	Multi wrench with ISO adapter	Non-sterile	Reusable
91116	Drill indicator for WS-75 and Osscora	Non-sterile	Reusable

Osscora Surgical Set Items Requiring Sterilization			
91102	Sterilisation cassette	Non-sterile	Reusable
91050	Osscora handpiece, 20:1	Non-sterile	Reusable
91054	Osscora motor with cable, 3.5m	Non-sterile	Reusable
04032600	Motor protective cover	Non-sterile	Reusable

Osscora main unit is also needed but not sterilized.

Required Sterile Items — New for Baha Attract

93550	BIM400 Implant magnet	Sterile	Single use
--------------	-----------------------	---------	------------

Required Sterile Items — Existing Baha Items Required for Baha Attract

92128	BI300 Implant, 3mm (for back-up)	Sterile	Single use
92129	BI300 Implant, 4mm	Sterile	Single use
93363	Conical guide drill, 3 and 4mm (92139 Guide drill, 3 and 4mm Ø 2.3mm if Conical Guide Drill not available)	Sterile	Single use
92140	Widening drill, 3mm with countersink (for back-up)	Sterile	Single use
92141	Widening drill, 4mm with countersink	Sterile	Single use
92136	Cover screw conical (for back-up)	Sterile	Single use

www.Cochlear.com/US

Cochlear Americas
13059 East Peakview Avenue
Centennial, CO 80111 USA
Telephone: 1 303 790 9010
Support: 1 800 483 3123

Cochlear Canada Inc.
2500-120 Adelaide Street West
Toronto, ON M5H 1T1 Canada
Telephone: 1 416 972 5082
Fax: 1 416 972 5083



©Cochlear Limited 2013. All rights reserved. Hear Now And Always and other trademarks and registered trademarks are the property of Cochlear Limited or Cochlear Bone Anchored Solutions AB. The names of actual companies and products mentioned herein may be the trademarks of their respective owners.

