

Cochlear™ Nucleus® CP802 Sound Processor

Product catalogue



Hear now. And always





INDEX	PAGE
CP802 Sound Processor	1
CP802 Battery options	4
CP802 Replacement parts	7
Retention options	8
Parental confidence	8
Care, maintenance and storage	9

Cochlear™ CP802 Sound Processor catalogue

In this catalogue you can find replacement parts and accessories for the CP802 Sound Processor. Your clinic can order the parts that you need. Simply fill out the order form and give it to your clinic.



COCHLEAR™ NUCLEUS®
CP802 SOUND PROCESSOR

Cochlear Nucleus CP802 Sound Processor

CP802 PROCESSING UNIT

The processing unit contains sophisticated micro-chip technology to process and deliver the sound to the implant and can be powered by various battery options.

Order number:
Z319055
(Carbon)



Z319056
(Maize)



CP802 COIL

The coil holds the coil magnet which sends the digital signals from the sound processor to the implant. Choose from two colours to complement your sound processor.

Order number:
Z324060
(Carbon)



Z324061
(Maize)




CP802 COIL CABLE

Coil cables connect the coil to the sound processor. Choose from two colours and four lengths to suit your individual needs.

Order number:

Z341661 (Carbon, 6cm)
Z341662 (Carbon, 8cm)
Z341663 (Carbon, 11cm)
Z341664 (Carbon, 28cm)



Z341665 (Maize, 6cm)
Z341666 (Maize, 8cm)
Z341667 (Maize, 11cm)
Z341668 (Maize, 28cm)



CP802 COIL MAGNET

Choose from two colours and six magnet strengths to suit your individual needs.

Order number:

Z341676 (Carbon, 1/2 M)
Z341677 (Carbon, 1 M)
Z341678 (Carbon, 2 M)
Z341679 (Carbon, 3 M)
Z341680 (Carbon, 4 M)
Z341681 (Carbon, 5 M)



Z341669 (Maize, 1/2 M)
Z341670 (Maize, 1 M)
Z341671 (Maize, 2 M)
Z341672 (Maize, 3 M)
Z341673 (Maize, 4 M)
Z341674 (Maize, 5 M)



Nucleus CP802 Battery options

CP802 STANDARD TAMPER RESISTANT BATTERY COVER

The disposable battery option offers long-lasting battery life of up to 60 hours and consists of a battery holder, tamper resistant cover and tamper resistant locking tool.

Order number:
Z319184
(Carbon)

Z319185
(Maize)



CP802 SOUND PROCESSOR BATTERY HOLDER


Order number:
Z198041



BATTERY COVER LOCKING TOOL

Secures the standard tamper resistant battery module to the processor to prevent your child from tampering with the unit.

Order number:
Z206768



ZINC AIR BATTERIES

Disposable batteries for use with the standard tamper resistant battery cover and battery holder.

Order number:
B13033



CP802 STANDARD RECHARGEABLE BATTERY MODULE

Rechargeable batteries are cost effective and good for the environment. The standard rechargeable battery module offers up to 31 hours of battery life. Choose from two colours

Order number:
Z319183
(Carbon)

Z319707
(Maize)



POWERBASE® BATTERY CHARGER C5015

The PowerBase charger is compact and includes four high quality rechargeable AAA batteries.

Order number:
Z60772



CP802 BATTERY PACK WITH CLIP

The battery pack fits into the battery clip retention case, allowing you to wear the batteries on the body, reducing the weight on the ears. The battery pack uses two AAA rechargeable or disposable batteries. Available in Carbon only.

Order number:
Z319186



CP802 BATTERY PACK CLIP

The battery pack clip is a retention case used to secure the Battery Pack to your clothes. Available in Carbon only.


Order number:
Z319180



CP802 BATTERY PACK CABLE

The battery pack cable connects the processing unit to the battery pack. Choose from two colours and three lengths to suit your individual needs.

Order number:
Z319088 (Carbon, 100cm)
Z319089 (Maize, 100cm)
Z319084 (Carbon, 30cm)
Z319085 (Maize, 30cm)
Z319086 (Carbon, 50cm)
Z319087 (Maize, 50cm)





CP800 SERIES SOUND PROCESSOR BATTERY CHARGER


Used for recharging battery modules
(recharge time takes 2-4 hours).

Order number:
Z341685



SOUND PROCESSOR CHARGING KIT

Order number:
Z342011 (Black*)



*Only available in black

Replacement parts

CP802 MICROPHONE COVERS

The microphone cover is designed to protect the microphones from dirt and moisture. Change the microphone cover if it looks dirty, or if you notice a gradual loss of sound quality. It is generally a good idea to change the microphone cover every two to three months.

Order number:
Z319699 (Carbon)

Z322676 (Maize)



CP802 MICROPHONE COVER SLEEVE

The microphone cover sleeve helps prevent children removing the microphone cover from the processing unit. The microphone cover sleeve is recommended for use by children under the age of 3 years.

Order number:
Z384719 (Carbon)

Z384720 (Maize)



CP802 TAMPER RESISTANT EARHOOK

This earhook comfortably and securely attaches the sound processor to your child's ear and can be moulded to suit different ear shapes and sizes. Tamper resistant locks ensure your child cannot remove the earhooks from the sound processor.

Order number:
Z319114 (Small)




CP802 EARHOOK

Earhooks comfortably secure the sound processor to your ear and can be moulded to ears of different shapes and sizes.

Order number:
Z319700 (Small)

Z319701
(Medium)

Z319702
(Large)




Retention options

MIC LOCK

The Mic Lock is a lightweight transparent plastic tube designed to firmly secure the BTE (behind the ear) processor to the ear. It can be cut to fit any size.

Order number:
Z300304



MIC LOCK STIRRUP

The mic lock stirrup is a lightweight transparent plastic tube designed to firmly secure the processor to the ear, while using the battery pack. It can be cut to fit any size.

Order number:
Z368868



Parental confidence

CP800 SERIES MONITOR EARPHONES

Ideal for when parents or teachers want to check the microphone input and sound quality on a child's sound processor.

Order number:
Z208300



Care, maintenance and storage

DRY & STORE®

This user friendly product removes moisture from your sound processor. It helps to increase the processor's longevity and may result in fewer repairs and better sound quality.

Order number:
Z60761



DRY-BRIK® DESICCANT

To be used in conjunction with the Dry & Store unit to remove moisture from your sound processor. Supplied as a pack of three.

Order number:
Z60764



STORAGE CASE

Store your processor, cables and accessories.

Order number:
Z342010





Hear now. And always

This is Cochlear's promise to you. As the global leader in hearing solutions, Cochlear is dedicated to bringing the gift of sound to people all over the world. With our hearing solutions, Cochlear has reconnected over 250,000 people to their families, friends and communities in more than 100 countries.

Along with the industry's largest investment in research and development, we continue to partner with leading international researchers and hearing professionals, ensuring that we are at the forefront in the science of hearing.

Our commitment to those who have received a Cochlear hearing solution is that for the rest of their life, they will Hear now. And always

Hear now. And always

Cochlear Ltd (ABN 96 002 618 073) 1 University Avenue, Macquarie University, NSW 2109, Australia Tel: 61 2 9428 6555 Fax: 61 2 9428 6352

Nihon Cochlear Co Ltd Ochanomizu-Motomachi Bldg, 2-3-7 Hongo, Bunkyo-Ku, Tokyo 113-0033, Japan Tel: 81 3 3817 0241 Fax: 81 3 3817 0245

Cochlear (HK) Ltd Unit 1810, Hopewell Centre, 183 Queens Road East, Wan Chai, Hong Kong SAR Tel: 852 2530 5773 Fax: 852 2530 5183

Cochlear Medical Device (Beijing) Co Ltd Unit 2208, Tower B, 91 Jianguo Road, Chaoyang District, Beijing 100022 P.R. China

Tel: 86 10 5909 7800 Fax: 86 10 5909 7900

Cochlear Limited (Singapore Branch) 6 Sin Ming Road, #01-16 Sin Ming Plaza Tower 2, Singapore 575585 Tel: 65 6553 3814 Fax: 65 6451 4105

Cochlear Korea Ltd 1st floor, Cheongwon building, 828-5, Yuksam dong, Kangnam gu, Seoul, Korea Tel: 82 2 533 4663 Fax: 82 2 533 8408

Cochlear Medical Device Company India Private Ltd Ground Floor, Platina Building, Plot No C-59, G-Block, Bandra Kurla Complex,

Bandra (E), Mumbai 400051 Tel: +91 22 61121111 Fax: +91 22 61121100

www.cochlear.com

The statements and images shown in this product catalogue are believed to be true and correct at the date of publication. However, specifications are subject to change without notice. Products shown in this catalogue may not be available in all countries. Please contact your local Cochlear office or Distributor for further information.

Cochlear, the elliptical logo, Nucleus and Hear now. And always are registered trademarks of Cochlear Limited.

Dry&Store is a registered trademark and Dry-Brik is a trademark of Ear Technology Corporation. © Cochlear Limited 2012.

N35823F ISS1 JUL07